Akku laden

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
Charging the battery

Fully charge your EndoLight® before first use: Connect the USB port of the charging cable to the power supply and plug it into a power outlet. The other end of the charging cable has a USB Type-C connector. Plug it into the socket on the EndoLight®.

Get started

The diodes contain sensors which are activated by skin contact: For irradiation of both arteries on your left wrist, ensure the diodes are resting on the skin. The EndoLight® must not be too tight, so as not to impede blood circulation (tight, but not constricting).

Warning

Never look directly into the diodes or shine the light into anyone else's eyes - visible and invisible radiation can be harmful to the eye.

Modes

Your EndoLight® has three modes that you can activate as needed. By repeatedly pressing the ... button, you can scroll through the modes. Press ⏻ to choose your selection.

Recharge

Activates both red and infrared diodes

Relax

Activates yellow, green and blue diodes

Recover

Activates all diodes simultaneously

EndoLight®

RECHARGE. RELAX. RECOVER.

by Weber Medical

Make your settings

⏻ Briefly press the On/Off button to switch on (a short press is sufficient)

⋮ Select one of the three modes

⏻ Confirm with On/Off

◀▶ Set intensity via rocker switch

⏻ Activate diodes via On/Off button

⏻ Switch off the device by pressing and holding the On/Off button or wait for timer to finish